

Title: Gilson Pipetman® Performance Verification SOP

Approvals:

Preparer: Judith Fitzpatrick Date 13Oct09
Reviewer: Kari Britt Date 13Oct09

1. Purpose:

To verify the calibration of a single channel pipette.

2. Scope:

Covers the cleaning, decontamination and verification of a single channel pipette.

3. Responsibilities:

3.1. It is the responsibility of the course instructor /lab assistant to ensure that this SOP is performed as directed and to update the procedure when necessary.

3.2. It is the responsibility of the students/technicians to follow the SOP as described and to inform the instructor about any deviations or problems that may occur while performing the procedure.

4. References:

4.1. Balance operation SOP

4.2. Balance calibration SOP

4.3. Autoclave SOP

4.4. Gilson Pipetman® Operation and Maintenance SOP

5. Definitions: N/A

6. Precautions: N/A

7. Materials:

7.1. Balance

7.2. 20g and 200g mass weights

7.3. Weigh boats

7.4. Deionized (DI) water

7.5. Small beaker for holding DI water

7.6. Verification labels

7.7. Verification form

7.8. Verification Pass/Fail form

7.9. Pipette tips

7.10. Pipetman® P20, P200, and P1000

7.11. 70% Isopropyl Alcohol (IPA)

7.12. Lab towels

7.13. Tweezers

7.14. Thermometer

7.15. Calculator

7.16. Barometer

Title: Gilson Pipetteman® Performance Verification SOP

8. Procedure:

8.1. Clean the pipette (See Figure 2 page 8)

Note: Most pipettes are designed so that the parts that normally come into contact with liquid contaminants can easily be cleaned and decontaminated.

- 8.1.1. Wipe entire pipette with a lab towel dampened with a mild detergent solution.
- 8.1.2. Wipe entire pipette with a lab towel dampened with distilled water.
- 8.1.3. Remove the tip ejector.
- 8.1.4. Wipe the tip ejector with a lab towel dampened with a mild soap solution.
- 8.1.5. Wipe the tip ejector with a lab towel dampened with distilled water.
- 8.1.6. Refit the tip ejector and allow the pipette to dry.
- 8.1.7. Dispose of lab towels in bio-hazardous waste receptacle.

8.2. Chemical decontamination

- 8.2.1. Spray a lab towel with 70% IPA to dampen the lab towel.
- 8.2.2. Wipe upper part of body with dampened lab towel.
- 8.2.3. Wipe tip holder tip ejector with dampened lab towel.
- 8.2.4. Wipe entire pipette with a lab towel dampened with distilled water.
- 8.2.5. Leave pipette to dry or wipe pipette dry with lab towel.
- 8.2.6. Dispose lab towels in bio-hazardous waste receptacle.

8.3. Verification of Calibration (See Figure 1 page 12)

Note: To test the accuracy of the pipette you will pipette a set volume 10 times and then weigh the total pipetted volume. 1mL of DI water should weigh 1g and 1µL should weigh 1mg. Calculate your % Error using the equation below:

$$\frac{\text{Expected Mass} - \text{Actual Mass}}{\text{Expected Mass}} \times 100 = \% \text{ Error}$$

If the % Error is $\leq 2\%$ the pipette passes verification if it is greater than 2% the pipette fails. We will verify the pipette once at the maximum volume for the pipette, once at the $\frac{1}{2}$ maximum volume, and once at the minimum volume. Altogether you will pipette 30 volumes and weigh 3 times for each pipette.

- 8.3.1. Record the necessary information on the Verification form.
- 8.3.2. Verify that the balance is still in calibration with a 20gm and 200g mass.
- 8.3.3. Verify that the calibration label of the balance is within the dated calibration time period.
- 8.3.4. Fill a small beaker with DI water.
- 8.3.5. Place the weigh boat on the balance.
- 8.3.6. Tare the balance and verify that 0.00 is being displayed.
- 8.3.7. Verify that the pipette is set to the maximum volume (e.g. the maximum volume for a P-20 pipette is 20µL.).

Title: Gilson Pipetteman® Performance Verification SOP

- 8.3.8. On the Pipette Verification Form, beside Selected Volume, enter the volume you will be pipetting and the value of that volume times 10 (e.g. for a 20 μ L pipette you will record 20 μ L for the selected volume and 200 μ L for the selected volume times 10.).
- 8.3.9. Calculate the expected mass by converting the selected volume times 10 using the following conversions: 1 μ L = 1mg and 1mL = 1g. Use the selected volume times 10 as the volume (e.g. for a 20 μ L pipette, 200 μ L multiplied by 1mg/ μ L = 200mg). Record the expected mass in the box beside Expected Mass.
- 8.3.10. Verify that the pipette is set to the maximum volume recommended by the manufacturer for the pipette.
- 8.3.11. Place pipette tip securely on the pipette.
- 8.3.12. Aspirate DI water into pipette tip from the beaker and dispense it into weigh boat. Refer to Figure 1 page 6 for proper operation of the pipette.
- 8.3.13. Repeat the above step 9 times. Each time you dispense the selected volume mark the Verification form in the numbered box beside Dispense Repetitions.
- 8.3.14. Record the final mass on the Verification form next to Recorded Mass.
- 8.3.15. Tare the balance and verify that 0.00 is being displayed.
- 8.3.16. Set the volume of the pipette to half capacity (e.g. For a P-20 pipette, set it to 10 μ L.) and verify the volume.
- 8.3.17. Repeat steps 8.3.9. through 8.3.16 with the pipette set to the half-capacity volume.
- 8.3.18. Tare the balance and verify that 0.00 is being displayed.
- 8.3.19. Set the volume of the pipette to the minimum capacity recommended by the manufacturer (e.g. For a P-20 pipette, set it to 2 μ L.)
- 8.3.20. Repeat steps 8.3.9. through 8.3.16 with the pipette set to the minimum-capacity volume.
- 8.3.21. Calculate the % Error (as directed in the note at the beginning of section 8.3) for each test (maximum, half-capacity, and minimum volumes) and record the results on the verification form.
- 8.3.22. Verify that all fields of the verification form have been filled out and fill out the Pipette Verification Pass/Fail form according to the results of the tests.

9. Attachments:

- 9.1. Figure 1: Calibration Equipment and Supplies
- 9.2. Pipette Performance Verification Form: Document 4.24.1
- 9.3. Pipette Performance Verification pass/Fail form: Document 6.24.2

10. History:

Name	Date	Amendment
Bob O'Brien	10Jul07	Initial Release
Bob O'Brien	11May08	College name change, add photographs, reword to clarify
Judith Fitzpatrick	13Oct09	Added directions for using 3 different volumes for calibrating the same pipette.

Title: Gilson Pipetteman® Performance Verification SOP

Mary Jane Kurtz and Sonia Wallman	28May10	Calibration changed to Performance Verification
--------------------------------------	---------	---



Di water
Waste Beaker
Pipette tips
Balance

Pipette rack
Tweezers
Small bottle

Figure 1: Calibration Equipment and Supplies