


**GFP Transformation Lab** 




#### **GFP Transformation Lab**









Images taken without permission from <a href="http://upload.wikimedia.org/wikipedihttp://bioinfo.biotec.or.th/Picture/Cell%20Tutorial/image005.jpg">http://bioinfo.biotec.or.th/Picture/Cell%20Tutorial/image005.jpg</a>, <a href="http://www.plantsci.cam.ac.uk/Haseloff/SITEGRAPHICS/Jellyfish.jpeg">http://www.plantsci.cam.ac.uk/Haseloff/SITEGRAPHICS/Jellyfish.jpeg</a>

## Bring Biotechnology to your Classroom

- Demonstrate the central framework of molecular biology
- Transform bacteria into glowing colonies
- Select for transformed Cells by antibiotic resistance
- View operon control over pGLO protein production
- Introduction to Biomanufacturing


## GFP (Green Fluorescent Protein

- Naturally produced in Jellyfish— Aequorea victoria
- Discovered in 1960's
- Source of bioluminescence when exposed to UV light



# Structure of the GFP Protein





Img Src: http://www.chem.leidenuniv.nl/metprot/armand/images/0291.jpg

### Detecting Gene Activity



- PGLO gene is inserted into DNA near a gene of interest
- It acts as a reporter gene
  - linked to another gene & glowing protein appears if it is expressed
- Expressed in entire animals







 $Img\ Src:\ http://www.bio.umass.edu/microscopy/images/gfp.jpg$ 



Img Src: http://www.mshri.on.ca/nagy/graphics/GFP%20mic e.jpg

Img Src: http://www.antville.org/img/pop/gfp.jpg


Img Src: http://www.computerra.ru/pubimages/73944.jpg



#### **Other Fluorescent Proteins**








## pGLO plasmid






- GFP gene
  - Codes for the GFP protein
- Bla gene
  - Codes for the enzyme  $\beta$ -lactamase
  - $\beta$ -lactamase destroys the antibiotic ampicillin
- araC regulator protein
  - Controls expression of GFP



#### Overall Goal of Lab Experiment




• Use genetic engineering techniques to insert the GFP gene into  $E.\ coli$ Plasmid containing gene of interest



## Selecting for Transformed Cells



- Selectable Marker: Trait that helps identify a transformed cell by conferring resistance to ampicillin
- Ampicillin presence in LBAgar will kill wild type E.coli BUT
- Transformed E. coli survive in the presence of ampicillin in LB Agar



Transformed E.coli


ampicillin

## **Arabinose Operon**

The arabinose operon in bacteria consists of the following.

Usually, the araC protein binds to the arabinose operon operator  $\rightarrow$  prevents transcription


When arabinose is present, it binds to the araC protein -> can't bind to operator → RNA polymerase can continue



Scientists modified the arabinose operon in pGLO to express the GFP gene.


araC protein binds to \_\_\_\_
the operator \(\frac{\rightarrow}{\rightarrow}\)
prevents transcription

When arabinose binds to araC it can no longer bind to operator → GFP gene is transcribed and translated





## Central Dogma of Molecular





- Spread E. coli without plasmid (- DNA) on plain LB agar
  - Wild type E. coli can grow demonstrated
- Spread E. coli without plasmid (- DNA) on LB/amp
  - E. coli aren't already resistant to ampicillin

### Transformation Yields Product



- What does this lead to?
  - Ability to produce a protein we need but can't make
  - Cell acts as the factory for the product under the correct conditions
  - Increased cell number yields increased product



### **Transformation Procedure**

- Step 1 Prepare appropriate plates
- Step 2 Suspend cells in CaCl<sub>2</sub> solution
- Step 3 Add pGLO plasmid to cells/put onto ice
- Step 4 Heat Shock at 42oC /put onto ice
- Step 5 Add nutrient broth to cells
- Step 6 Streak cells on to appropriate plates



#### Transformation Time Line


- First step: Grow up colonies of E.coli
- Second step:
   Prepare Selective
   media
- Transform cells with pGLO plasmid
- Detect transformed cells

- 2-3 days required
- 1 day
- 45 minutes
- Results in 24 hours
- Supplies for up to 32 students

### PGLO Transformed E.coli



- Cells containing pGLO plasmid are now resistant to ampicillin
- Cells containing pGLO plasmid will also glow green when arabinose



#### Biomanufacturing



 Upstream Processing: Growing genetically transformed cells that produce a desired protein

 Downstream Processing: Separation and purification of that product for human use