

SOP: Cryopreservation of NISTCHO Cells

Approvals:

Preparer: Hetal Doshi

Date: 27SEP2024

Reviewer: Dr. Maggie Bryans

Date: 01NOV2024

1. Purpose:

- 1.1. This Standard Operating Procedure (SOP) describes the steps required for the cryopreservation of NISTCHO cells. NISTCHO cells are a Chinese hamster ovary (CHO) cell line that has been engineered by MilliporeSigma to express a recombinant anti RSV monoclonal antibody (cNISTmAb)

2. Scope:

- 2.1. This SOP will be applied to cryopreservation of NISTCHO in suspension when cell culture reaches a density of $3 \times 10^6 - 4 \times 10^6$ viable cells/mL (SOP: Resuscitation and Culture of NISTCHO cells)

3. Summary of Method:

- 3.1. Determine the cell concentration and viability of the NISTCHO culture in the log phase.
- 3.2. Calculate and prepare the amount of cryopreservation medium required to cryopreserve the cells at 1.3×10^7 cells/ml concentration.
- 3.3. Centrifuge the appropriate amount of cell suspension and resuspend them in a appropriate volume of cryopreservation medium
- 3.4. Place the vials in freezing cube/Mr. frosty and store at -80°C overnight and transfer to the -150°C freezer or LN2 freezer.

4. References:

- 4.1. NISTCHO Test Material, Clonal CHO-K1 Cell Line Producing cNISTmAb Guidance Document <https://tsapps.nist.gov/srmext/certificates/10197.pdf>
- 4.2. SOP: Labconco Purifier Class 2 Biological Safety Cabinet Operation, Document No. UP 1
- 4.3. SOP: Operation of Logos biosystems Luna-FL Fluorescence Cell Counter for Fluorescence Cell Counting, Document No. UP22
- 4.4. SOP: Resuscitation and Culture of NISTCHO Cells

5. Definitions:

- 5.1. N/A

6. Precautions:

- 6.1. Use BSL1 safety measures and practices and discard waste in biohazard containers
- 6.2. Routine care should be exercised in the handling of buffers and samples of biological materials, which may have harmful biological activity in the case of accidental ingestion, needle stick etc.
- 6.3. Gloves, a lab coat and protective eyewear should be worn

7. Responsibilities:

- 7.1. It is the responsibility of the course instructor/lab assistant to ensure that this SOP is performed as described and to update the procedure when necessary.
- 7.2. It is the responsibility of the students/technician to follow the SOP as described and to inform the instructor about any deviations or problems that may occur while performing the procedure.

8. Equipment and Materials:

8.1. Equipment:

SOP: Cryopreservation of NISTCHO Cells

- 8.1.1. Biological safety cabinet Class II (BSC)
- 8.1.2. CO₂ incubator with shaking capacity set at 37°C ±2, 5% CO₂, and 125rpm (e.g., Benchmark incu-shaker™ CO2 mini, with non-slip rubber mat platform catalog# sku: H3501)
- 8.1.3. Eppendorf centrifuge 5804R
- 8.1.4. Luna-FL Fluorescence Cell Counter
- 8.1.5. Vacuum pump
- 8.1.6. Ultra-Low temperature freezer
- 8.1.7. Controlled Freeze Cubes, Fisherbrand™ Catalog # 13-131-014

8.2. Materials:

- 8.2.1. NISTCHO cell culture 3 X 10⁶ cells/ml- 4 X 10⁶ cells/ml cell concentration with a viability of ≥ 95%
- 8.2.2. EX-CELL® CD CHO Fusion medium, Sigma Aldrich catalog #14365C-1000ML, (Will be referred as expansion medium)
- 8.2.3. Dimethyl Sulfoxide, Fisher Bioreagents, Catalog# BP231-100
- 8.2.4. Cryovials, Corning, Catalog # 8670
- 8.2.5. Pipette aid.
- 8.2.6. Cryovial rack
- 8.2.7. Sterile serological pipettes (2ml, 5ml and 25 ml)
- 8.2.8. Lab coat, gloves, sleeves
- 8.2.9. 70% Isopropanol or Ethanol
- 8.2.10. Lint free clean room wipes
- 8.2.11. Kim Wipes
- 8.2.12. 1.5 ml microfuge tube and tube holder
- 8.2.13. 15ml conical tube with stand
- 8.2.14. 15ml conical tube with appropriate volume of milliQ water to use as a blank for centrifuge.
- 8.2.15. Autoclaved glass Pasteur pipette
- 8.2.16. 125ml shake flask holder (magnetic)
- 8.2.17. Acridine Orange/Propidium Iodide Stain Logos Biosystem, Catalog # F23001
- 8.2.18. Cryovial marker
- 8.2.19. 0.22-micron PES membrane syringe filter
- 8.2.20. Sterile syringe

9. Procedure:

- 9.1. Determination of cell concentration and cell viability of the existing NISTCHO cell culture
 - 9.1.1. Prepare the BSC as per the SOP and verify the controlled freezer cube is at room temperature.
 - 9.1.2. Collect the following material and place it in the biological safety cabinet after swabbing or wiping them with 70% Ethanol.
 - Sterile serological pipette size 1ml
 - Pipette aid
 - 1 1.5 ml sterile microfuge tube

SOP: Cryopreservation of NISTCHO Cells

- Cryovial rack
- 9.1.3. Remove shake flask with NISTCHO cells from 37°C incubator with 5% CO₂ and place the flask in the BSC after swabbing it with 70% Isopropanol. Aseptically transfer 200 µl of the cell suspension to a sterile microfuge tube for determination of viable cell count and cell viability.
- 9.1.4. Place the shake flask back into the 37°C incubator with 5% CO₂ shaking at 125 rpm.
- 9.1.5. Determine the viable cell density and viability using Luna-FL Fluorescence Cell Counter
- 9.2. Preparation of Cryopreservation Medium
- 9.2.1. Calculate the volume of cryopreservation medium needed for cryopreservation of cells to a final density of 1.3×10^7 cells/ml/vial (Volume based on total number of cells available for freezing)
- 9.2.2. Gather the following items, spray or wipe with 70% Isopropanol, and place in the biological safety cabinet.
- Pack of sterile 15ml conical tube
 - Conical tube rack/holders
 - Cryovial rack
 - Sterile serological pipette sizes 2ml, 5ml, 10ml, and 25ml
 - Pipette aid.
 - Expansion Medium (EX-CELL® CD CHO Fusion medium, Sigma Aldrich catalog #14365C-1000ML)
 - Dimethyl Sulfoxide
 - Autoclaved glass Pasteur pipette (if not already in the BSC)
 - Sterile syringe
 - 0.22-micron Syringe filter
 - Pack of sterile 50 ml conical tube
 - Conical tube with appropriate amount of milliQ water to be used as a blank for centrifuge
- 9.2.3. Prepare the appropriate volume of cryopreservation medium with 93% Expansion Medium+ 7% DMSO. (volume calculated in step 9.2.1. plus 2ml)
- 9.2.4. Sterile filter the prepared cryopreservation medium with 0.2 µm syringe filter or bottle top filter depending on volume
- 9.2.5. In BSC label appropriate number of vials as “NISTCHO cells, Passage number, cell concentration/ml/vial, Date, and initials”
- 9.2.6. Remove shake flask with NISTCHO cells from 37°C incubator with 5% CO₂ and place the flask in the BSC after swabbing it with 70% Isopropanol.
- 9.2.7. Aseptically transfer the appropriate volume of cell suspension to the conical tube
- 9.2.8. Balance tube and centrifuge at 220 x g for 5 minutes at 20°C
- 9.2.9. Place the tube in BSC after spraying with 70% Isopropanol. In the BSC discard the supernatant without disturbing cell pellet
- 9.2.10. Resuspend pellet in 10% of final banking volume of cryopreservation medium prepared in step 9.1.5.

SOP: Cryopreservation of NISTCHO Cells

- 9.2.11. Adjust the volume with cryopreservation medium to obtain the desired cell density of 1.3×10^7 cells/ml
- 9.2.12. Aliquot 1mL of cell suspension into 2mL screw-capped cryovials. Apply the cap to the cryovials, seal well.
- 9.2.13. Place the vials in the Controlled Freeze Cube
- 9.3. Freeze the cells at -80°C overnight prior to transferring the frozen cells to -150°C .

10. History:

Revision Number	Effective Date	Preparer	Description of Change
0	01NOV2024	Hetal Doshi	Initial release