

SOP: Cryopreservation of Suspension Viral Production Cells 2.0 - HEK293F

Approvals:

Preparer: Isso Bayala & Hetal Doshi
Reviewer: Dr. Maggie Bryans

Date: 15AUG22
Date: 16AUG22

1. Purpose:

- 1.1. This Standard Operating Procedure (SOP) describes the steps required for the cryopreservation of Viral Production Cells 2.0 in suspension under BSL-2 safety criteria. Viral production Cells 2.0 Prototype are a clonal derivative of HEK293F cell line and have been adapted to suspension, high-density culture in Gibco Viral Production Medium. These cells can be thawed directly into Gibco Viral Production Medium

2. Scope and Applicability:

- 2.1. This SOP will be applied to cryopreservation of viral production cells 2.0 cells in suspension when cell culture reaches a density of 4×10^6 - 6×10^6 viable cells/mL (SOP: Resuscitation and Culture of Viral Production Cells (VPC) 2.0 -HEK293F UP 32)

3. Responsibilities

- 3.1. It is the responsibility of the course instructor/lab assistant to ensure that this SOP is performed as described and to update the procedure when necessary.
- 3.2. It is the responsibility of the students/technician to follow the SOP as described and to inform the instructor about any deviations or problems that may occur while performing the procedure.

4. References:

- 4.1. Clonal HEK293F-derived cell line adapted for AAV production user manual
https://www.thermofisher.com/document-connect/document-connect.html?url=https%3A%2F%2Fassets.thermofisher.com%2FTFS-Assets%2FSLG%2Fmanuals%2FMAN0019620_ViralProductionCells_2-0_and_ViralProductionMedium_UG.pdf
- 4.2. SOP: Labconco Purifier Class II Biological Safety Cabinet Operation, Document No. UP 1
- 4.3. SOP: Bio-Rad TC20 Cell Counter
- 4.4. SOP: Trypan Blue Assay, Document No. UP6

5. Precautions

- 5.1. Use BL2 safety measures and discard waste in biohazard containers.
- 5.2. Routine care should be exercised in the handling of buffers and samples of biological materials, which may have harmful biological activity in the case of accidental ingestion, needle stick etc
- 5.3. Gloves, a lab coat and protective eyewear should be worn when handling buffers and samples.

6. Equipment and Materials:

- 6.1. Equipment
 - 6.1.1. Biological safety cabinet
 - 6.1.2. Liquid Nitrogen Dewar
 - 6.1.3. Equipment to determine cell viability (cell counter or hemocytometer)
 - 6.1.4. centrifuge
 - 6.1.5. Compound Light Microscope with 100X magnification (10X objective lens)
 - 6.1.6. Ultra-Low temperature freezer

SOP: Cryopreservation of Suspension Viral Production Cells 2.0 - HEK293F

6.2. Materials

- 6.2.1. Viral Production Cell (VPC) Culture
- 6.2.2. 2.0 mL screw-capped cryovials.
- 6.2.3. Gibco viral production Medium A49842DK
- 6.2.4. Gibco GlutaMAX Supplement (100X)35050061
- 6.2.5. HPLC grade dimethyl sulfoxide (DMSO)
- 6.2.6. Sterile 50-mL centrifuge tubes
- 6.2.7. Pipette aid
- 6.2.8. Sterile serological pipettes (2ml, 5ml and 25ml)
- 6.2.9. Lab coat, gloves, sleeves
- 6.2.10. 70% Isopropanol
- 6.2.11. Sterile cleaning wipes
- 6.2.12. Trypan Blue (0.4% solution)
- 6.2.13. 1.5 ml microfuge tube and tube holder
- 6.2.14. P20 micropipettes and compatible tips

7. Procedure:

7.1. Preparation of freezing Media

- 7.1.1. Prepare biological safety cabinet per Labconco Purifier Class 2 Biological Safety cabinet (BSC) Operation SOP
- 7.1.2. Gather the following items, spray or wipe with 70% Isopropanol, and place in the biological safety cabinet.
 - Pipette aid (sanitize with cleaning wipes or 70% IPA)
 - 25ml sterile pipettes
 - 5ml sterile pipettes
 - 1 ml sterile pipette
 - 125-mL PETG Erlenmeyer shaker flask
 - Cryovial rack
 - Gibco viral production Medium
 - GlutaMax supplement 200mM
 - DMSO
 - 1.5 ml Eppendorf tube
- 7.1.3. Prepare the freezing medium with 90% Viral Production Medium supplemented with 4mM GlutaMAX™ Supplement + 10% DMSO.
- 7.1.4. Sterile filter the prepared media with 0.2 µm syringe filter and keep on ice until ready to use

- 7.2. Remove VPCs 2.0 shaker flask from 37°C incubator with 8% CO₂ and place the flask in the BSC after swabbing it with 70% Isopropanol. Aseptically transfer remove 200 µl of the cell suspension for determination of viable cell count and cell viability
- 7.3. Place the shaker flask back into the 37°C incubator with 8% CO₂
- 7.4. Determine the viable cell density and viability using automated cell counter

SOP: Cryopreservation of Suspension Viral Production Cells 2.0 - HEK293F

- 7.5. Calculate the volume of freezing media needed for cryopreservation of cells to a final density of 1×10^7 cells/ml
- 7.6. Remove VPCs 2.0 shaker flask from 37°C incubator with 8% CO₂ and place the flask in the BSC after swabbing with 70% Isopropanol. Aseptically transfer appropriate amount of cell suspension into a sterile 50mL centrifuge tube in the BSC
- 7.7. Balance tube and centrifuge at $300 \times g$ for 5 minutes.
- 7.8. Place the tube in BSC after spraying with 70% Isopropanol. In the BSC discard the supernatant without disturbing cell pellet.
- 7.9. Resuspend pellet in 10% of final banking volume of cryopreservation medium prepared in step 7.1.3.
- 7.10. Adjust the volume with cryopreservation medium to obtain the desired cell density of 1×10^7 cells/ml.
- 7.11. Aliquot 1-mL of cell suspension into 2.0 mL screw-capped cryovials. Apply the cap to the cryovials, seal well
- 7.12. Place the vials in Mr. Frosty
- 7.13. Freeze the cells at -80°C for overnight prior to transfer the frozen cells to liquid nitrogen.
- 7.14. Transfer frozen the vials to liquid nitrogen for long-term storage

8. History:

Revision Number	Effective Date	Preparer	Description of Change
0	13APR22	Isso Bayala, Hetal Doshi	Initial release