

SOP: SDS-PAGE

Approvals

Preparer:	Lara Dowland	Date:	19Jan10
Reviewer:	Cianna Cooper	Date:	13Jun18

1. Purpose:

- 1.1. To describe the appropriate operating instructions to perform SDS PAGE analysis of proteins samples.

2. Scope:

- 2.1. Applies to confirming the presence and purity of protein samples.

3. Responsibilities:

- 3.1. It is the responsibility of the course instructor/lab assistant to ensure that this SOP is performed as described and to update the procedure when necessary.
- 3.2. It is the responsibility of the students/technicians to follow the SOP as described and to inform the instructor about any deviations or problems that may occur while performing the procedure.

4. References:

- 4.1. Invitrogen Novex Gel instructions
- 4.2. Fisher Scientific Vertical MiniGel System SOP

5. Definitions: N/A

6. Precautions:

- 6.1. Acrylamide is a neurotoxin. Always wear protective gloves when handling the polyacrylamide gels.

7. Materials:

- 7.1. Protein Samples
- 7.2. Precision Plus Protein Standards, Kaleidoscope from BioRad (2– 8°C refrigerator)
- 7.3. Fisher Scientific Vertical MiniGel System
- 7.4. Power Supply for Protein Electrophoresis
- 7.5. NuPAGE 4-12% Bis-Tris Gels (1.0mm x 10 well), 2– 8°C refrigerator
- 7.6. NuPAGE MOPS SDS Running Buffer (20X), room temperature
- 7.7. NuPAGE Antioxidant, 2– 8°C refrigerator
- 7.8. NuPAGE LDS Sample Buffer (4X), room temperature
- 7.9. Reducing Agent (10X), -20°C freezer
- 7.10. Graduated cylinders (1L and 100ml)
- 7.11. P20, P100 or P200 Micropipettor and tips, including gel loading tips
- 7.12. Microfuge Tubes
- 7.13. Microfuge
- 7.14. Heat block
- 7.15. Staining Trays
- 7.16. rotary shaker
- 7.17. Invitrogen Simply Blue Safe Stain
- 7.18. Light Box
- 7.19. Gel Documentation Instrument

8. Procedure:

SOP: SDS-PAGE

8.1. Prepare Running Buffer and Staining Solutions

- 8.1.1. 1L NuPAGE MOPS SDS Running Buffer (1X) (if needed)
 - 8.1.1.1. Place 50ml NuPAGE MOPS SDS Running Buffer (20X) in a 1 Liter graduated cylinder.
 - 8.1.1.2. Gently add 950ml deionized water by running it down the side of the cylinder to make 1 liter of 1X NuPAGE MOPS SDS Running Buffer.
 - 8.1.1.3. Add a stir bar and gently stir. NOTE: SDS is a detergent and will foam if mixed vigorously.
- 8.1.2. 200ml NuPAGE MOPS SDS Running Buffer (1X) plus antioxidant (if needed)
 - 8.1.2.1. Add 200ml of 1X NuPAGE MOPS SDS Running Buffer into a 500ml Erlenmeyer flask.
 - 8.1.2.2. Add 500µl of NuPAGE Antioxidant.
 - 8.1.2.3. Add a stir bar and gently stir. NOTE: SDS is a detergent and will foam if mixed vigorously.
 - 8.1.2.4. Store at 2– 8°C

8.2. Prepare Protein Samples

(Do NOT perform this step with the Molecular Weight Marker)

- 8.2.1. Dilute conditioned media from cell culture or flow through fraction off the chromatography column 1:4 in 1x PBS as follows:
 - 8.2.1.1. Pipette 9 ul of 1x PBS into 1.5ml centrifuge tube.
 - 8.2.1.2. Add 3 ul of sample to the same tube.
- 8.2.2. For **all** the samples, combine the following in a sterile 1.5ml microfuge tube:
 - 2.5ul 4x sample buffer
 - 1.0ul 10x reducing agent
 - 5ul sample
 - 1.5ul deionized water
- 8.2.3. Mix gently with a pipet by aspirating and dispensing at least 3 times
- 8.2.4. Heat at 70°C for 10 minutes.
- 8.2.5. Pulse all samples and standards in a microfuge for 30 seconds.

8.3. Prepare Fisher Scientific Vertical MiniGel Gel Box

- 8.3.1. Assemble gel box according to the Fisher Scientific Vertical MiniGel System SOP.
- 8.3.2. Place 200ml NuPAGE MOPS SDS Running Buffer (1X) plus antioxidant in upper buffer chamber (small chamber between 2 gels or the gel and buffer dam)
- 8.3.3. Fill the lower buffer chamber with approximately 600ml of 1X NuPAGE MOPS SDS Running Buffer (large chamber).
- 8.3.4. Rinse gel wells with micropipettor and buffer from upper buffer chamber.

8.4. Load Samples

- 8.4.1. Using a micropipettor and gel loading tips, load 10ul of the Molecular Weight Marker into one well and 10µl of each sample into separate wells.
 - 8.4.1.1. Avoid loading samples symmetrically.
- 8.4.2. Load any empty wells with 15µl of diluted 4X LDS Sample Buffer.
- 8.4.3. Record order of samples and volumes loaded.

SOP: SDS-PAGE

8.5. Run NOVEX NuPAGE MOPS SDS Precast Gel

- 8.5.1. Plug electrophoresis chamber into the gel electrophoresis power supply.
- 8.5.2. Run gel at 125V for 1 hour.
- 8.5.3. Turn off the power supply when the dye reaches 1cm from the bottom of the gel.

8.6. Stain and Photodocument the NOVEX NuPAGE MOPS SDS Precast Gel

- 8.6.1. Disassemble gel box per Fisher Scientific Vertical MiniGel System SOP and remove gel from plastic cassette.
- 8.6.2. Rinse gel box well with DI water. Do not use brushes on the gel box, they scratch the surface. Do not immerse top of gel box or electrical components.
- 8.6.3. Place gel in staining tray with 100ml of deionized water for 5 minutes with gentle shaking.
- 8.6.4. Repeat wash step 2x for a total of 3 washes.
- 8.6.5. Add 20ml of SimplyBlue SafeStain, cover with saran wrap, and gently shake for 1 hour at room temperature.
- 8.6.6. Wash in 100ml of deionized water for 1 hour or overnight with gentle shaking.
- 8.6.7. Remove gel from staining tray and place on visible light box
- 8.6.8. Identify the protein standards and samples and estimate their molecular weights.
See Molecular Weight Diagram.

9. Attachments:

- 9.1. Molecular Weight Marker Diagram

10. History:

Name	Date	Amendment
Lara Dowland	19Jan10	Initial Release
Lara Dowland	13Jun18	Update running conditions and staining procedure.

SOP: SDS-PAGE

