

SOP: CHO DP12 Spinner Flask Harvest

Approvals:

Preparer: Robin Zuck
Reviewer: Hetal Doshi
Reviewer: Dr. Maggie Bryans

Date: 25MAY2021
Date: 27MAY2021
Date: 02JUN2021

1. Purpose:

- 1.1. To harvest anti IL-8 mAb containing conditioned media from a Spinner Flask culture and prepare the harvested media for protein A chromatography using a gravity flow column.

2. Scope:

- 2.1. This SOP applies to the preparation of CHO DP12 conditioned media for the isolation of produced anti IL-8 monoclonal antibody using a gravity flow protein A column.

3. Summary of Method:

- 3.1. Transfer of the cell suspension from the spinner flask to 50ml centrifuge tubes.
- 3.2. Pellet the cells by centrifugation.
- 3.3. Clarify media by Microfiltration

4. References:

5. Responsibilities:

- 5.1. It is the responsibility of the course instructor/lab assistant to ensure that this SOP is performed as described and to update the procedure when necessary.
- 5.2. It is the responsibility of the students/technician to follow the SOP as described and to inform the instructor about any deviations or problems that may occur while performing the procedure.

6. Equipment and Materials:

- 6.1. 100 ml Spinner flask culture cell 6×10^5 - 1×10^6 cells/ml or 2 T75 flasks 80-100% confluent
- 6.2. 50 mL centrifuge tubes (2)
- 6.3. 0.2 micron vacuum sterile filtration unit with storage bottle, 100ml or 250ml bottle.
- 6.4. 25ml or 50ml serological pipettes and pipette aid
- 6.5. 70% Ethanol solution
- 6.6. Laboratory wipes
- 6.7. 37°C, 5% CO₂ Incubator
- 6.8. Biological Safety Cabinet
- 6.9. Benchtop Centrifuge

7. Procedure:

- 7.1. Transfer of culture from spinner flask to centrifuge tubes.
 - 7.1.1. Remove the Spinner Flask or T75 flasks from the incubator and place in the BioSafety cabinet.

SOP: CHO DP12 Spinner Flask Harvest

- 7.1.2. Aseptically transfer the culture from the Spinner Flask or T75 flasks to one or two - 50 mL centrifuge tubes using 25ml or 50 mL serological pipets and a PipetAid.
- 7.2. Centrifuge the culture in a pre-chilled centrifuge and rotor at 1000 x g for 5 minutes, at 4 ° C.
- 7.3. Place the centrifuge tubes in the Biosafety cabinet after swabbing with 70% ethanol. Aseptically transfer the supernatant, the conditioned medium (CM), from the centrifuge tubes to a 100ml or 250ml 0.22µm sterile filter unit with storage bottle by carefully decanting the supernatant off the cell pellet or pipetting the supernatant using a 25ml or 50ml pipette.
- 7.4. Filter the supernatant by applying a vacuum. Remove the filter portion of the unit and cap the bottle. Label the bottle “Anti-IL8 mAb conditioned media ,[Date], [final cell concentration],[Initial], [volume].Store the Conditioned Media at 4°C for further processing

9. Attachments:

10. History

<i>Revision Number</i>	<i>Effective Date</i>	<i>Preparer</i>	<i>Description of Change</i>
0		Robin Zuck	Initial release